

การพัฒนาแกรนูลชีไอพีลสีน่องกุทธินานศ์วัย เอเชียติกเรเดต เมซิลเมชาคริเตต
ไคไฟลิเมอร์ และเอเชียติกเรดจูไ/do:ชนิดกระจาดตัวในน้ำ^๑
โดยเทคนิคฟลูอิดไดซ์ เบค

นางสาวศรภาพ หริการภักดี



สถาบันวิทยบริการ
จุฬาลงกรณ์มหาวิทยาลัย
วิทยานิพนธ์นี้เป็นส่วนหนึ่งของการศึกษาตามหลักสูตรปริญญาเกสัชศาสตรมหาบัณฑิต
ภาควิชาเกสัชฯ สาขาวรรณ
บัณฑิตวิทยาลัย จุฬาลงกรณ์มหาวิทยาลัย
ปีการศึกษา 2539

ISBN 974-634-999-6

ลิขสิทธิ์ของบัณฑิตวิทยาลัย จุฬาลงกรณ์มหาวิทยาลัย

DEVELOPMENT OF THEOPHYLLINE SUSTAINED RELEASE GRANULES
COATED BY FLUIDIZED BED TECHNIQUE WITH POLY (ETHYLACRYLATE
METHYLMETHACRYLATE) AND ETHYLCELLULOSE
AQUEOUS DISPERSIONS

MISS SARAPORN HARIKARNPAKDEE

สถาบันวิทยบริการ
จุฬาลงกรณ์มหาวิทยาลัย
A Thesis Submitted in Partial Fulfillment of the
Requirements for the Degree of Master of Science in Pharmacy

Department of Manufacturing Pharmacy

Graduate School

Chulalongkorn University

Academic Year 1996

ISBN 974-634-999-6

Thesis Title Development of Theophylline Sustained Release Granules
Coated by Fluidized Bed Technique with Poly (ethylacrylate
Methylmethacrylate) and Ethylcellulose Aqueous Dispersions

By Miss Saraporn Harikarnpakdee

Department Manufacturing Pharmacy

Thesis Advisor Associate Professor Garnpimol C. Ritthidej, Ph.D.

Thesis Co-Advisor Assistant Professor Wichein Thanindrata, M.Sc. in Pharm

Accepted by the Graduate School, Chulalongkorn University in Partial
Fulfillment of the Requirements for the Master's Degree

Santi Thungsawan Dean of Graduate School
(Associate Professor Santi Thungsawan, Ph.D.)

Thesis Committee


P. Kulvanich Chairman
(Associate Professor Poj Kulvanich, Ph.D.)

Garnpimol C. Ritthidej Thesis Advisor
(Associate Professor Garnpimol C. Ritthidej, Ph.D.)

Wichein Thanindrata Thesis Co-Advisor
(Assistant Professor Wichein Thanindrata, M.Sc. in Pharm)

Ubonthip Nimmannit Member
(Associate Professor Ubonthip Nimmannit, Ph.D.)

Sirisak Dumrongpisudthigul Member
(Assistant Professor Sirisak Dumrongpisudthigul, M.Sc. in Pharm)



คราวพร หรือการภักดิ : การพัฒนาแกรนูลชิโอฟลิดลินออกฤทธิ์นานคัวว์ เอชิโอะคริเกต เมซิດ เมชาคริเกต ໄโค โพลิเมอร์ และเอชิಡเซกู ໄอกชนิดกระษายตัวในน้ำ ໄดบเทกนิกฟลูอิດไบซ์ เบค (DEVELOPMENT OF THEOPHYLLINE SUSTAINED RELEASE GRANULES COATED BY FLUIDIZED BED TECHNIQUE WITH POLY(ETHYLACRYLATE METHYLMETHACRYLATE) AND ETHYLCELLULOSE AQUEOUS DISPERSIONS) อ.ที่ปรึกษา : รศ.ดร.กาญจน์พิม พุทธิเดช, อ.ที่ปรึกษาร่วม : ผศ.ดร. วิเชียร ชานินทร์ราชาร, 232 หน้า. ISBN 974-634-999-6

การเตรียมชิโอฟลิดลินแกรนูลขนาดต่างๆ เกิดขึ้นด้วยเอชิಡเซกู ໄอกชนิด เอเชวิส ดีสเพอชัน ชุ่มน้ำแกป์ ไอซิตาส์ โดยใช้เทคนิคการหันเคลื่อนแบบฟลูอิດไบซ์ เบค เพื่อศึกษาคุณสมบัติทางกายภาพ และ ฤทธิ์สมบัติการปลดปล่อยตัวจากแกรนูลที่ถูกเกลื่อนแล้วซึ่งบรรจุในแคปซูล เปรียบเทียบกับแกรนูลที่ไม่ได้เกลื่อน พื้นผิวดองแกรนูลที่ไม่ได้เกลื่อนจะมีลักษณะหนาแน่นมากกว่าและมีรูปร่างไม่สม่ำเสมอ ส่วนพื้นผิวดองแกรนูลที่เกลื่อนด้วยบุคคลิค เอ็นอี 30°C จะมีลักษณะหนาแน่นกว่าพื้นผิวดองแกรนูลที่เกลื่อนด้วยชัวร์.ลีส การเพิ่มระดับการเกลื่อนจะทำให้ความหนาของพื้นผิวน้ำหนาขึ้น ดังนั้นแกรนูลจึงมีขนาดใหญ่ขึ้น รวมทั้ง ขนาดและนุ่มนวลต่อการดัดแปลงด้วยเมื่อมีการเกลื่อน ค่าแรงกิล ของทริพส แต่ค่าปริมาณความชื้นของ แกรนูลทั้งหมดจะอยู่ในช่วง 30-38 องศา และ 0.74-0.86 เปอร์เซ็นต์ ตามลำดับ การเปรียบเทียบ เปอร์เซ็นต์ก้อนเพรสซิบลิตี้ และค่าอัตราการไขтратะว่างแกรนูลที่ไม่ได้เกลื่อนและแกรนูลที่เกลื่อน แล้วในแต่ละขนาดของแกรนูล จะให้ผลการทดสอบที่มีค่าแปรปรวน เปอร์เซ็นต์อิส丢了เกชั่น แยก เปราก และเปอร์เซ็นต์ต่อเทอร์ ของปริมาณชัวร์ลีสพิล์นจะมีค่าเพิ่มสูงขึ้น เมื่อเปอร์เซ็นต์ของไคลบิวชัน พากาเดตที่ใช้เป็นพลาสติกใช้เซอร์วิคั่นสูงขึ้น ในขณะที่ค่าเทียนชาด์ สเกริงของพื้นจะลดต่อไป แต่ค่าปริมาณของเอเชวิส โพลิเมอร์ลิก ดีสเพอชัน ชนิดของเอเชวิสโพลิเมอร์ลิก ดีสเพอชัน จะมีผลต่อ การปลดปล่อยตัวยา แกรนูลที่มีขนาดใหญ่กว่าซึ่งบรรจุในแคปซูลจะให้การปลดปล่อยตัวยาที่ช้ากว่า แกรนูลที่มีขนาดเล็กกว่า ในปริมาณการปลดปล่อยตัวยาที่ใกล้เคียงกัน จะใช้ปริมาณของบุคคลิค เอ็นอี 20°C ที่ใช้ในการเกลื่อนสูงกว่าปริมาณของชัวร์ลีส รูปแบบการปลดปล่อยตัวจากแกรนูลที่บรรจุใน แคปซูลจะเป็นแบบเพร์สอยเดอร์และแกรนูลที่เกลื่อนแล้วบางครั้งทำรับจะให้การปลดปล่อยตัวยาที่เข้า นาตราเจนของยูอีสพี

C775247 : MAJOR MANUFACTURING PHARMACY

KEY WORD: ETHYLCELLULOSE / POLY(ETHYLACRYLATE METHYLMETHACRYLATE)

SARAPORN HARIKARNPAKDEE : DEVELOPMENT OF THEOPHYLLINE SUSTAINED RELEASE GRANULES COATED BY FLUIDIZED BED TECHNIQUE WITH POLY(ETHYLACRYLATE METHYLMETHACRYLATE) AND ETHYLCELLULOSE AQUEOUSE DISPERSION. THESIS ADVISOR : ASSO.PROF.GARNPIMOL C. RITTHIDEJ, Ph.D., THESIS CO-ADVISOR : ASSIS. PROF. WICHEIN THANINDRATARN, 232 pp. ISBN 974-634-999-6

Theophylline granules of different sizes were coated with ethylcellulose aqueous dispersion (Surelease[®]) or poly(ethylacrylate methylmethacrylate) aqueous dispersion (Eudragit[®] NE 30D) containing cab-o-sil by fluidized bed technique. The physical properties of coated granules and sustained release characteristics from capsules containing coated granules were evaluated and compared to those of uncoated granules. All uncoated theophylline granules exhibited rough surface and irregular shape. The granules coated with Eudragit[®] NE 30D exhibited rougher surface than Surelease[®] coated granules. Increasing the coating level increased the thickness of film, thus granules became slightly larger. Edge and corner were lower when coated.

The angles of repose and the moisture contents of all granules were within the range of 30-38° and 0.74-0.86 % respectively. Comparison of the percent compressibility and the flow rate between uncoated and coated granules of each size showed varied results. For Surelease[®] film, the percent elongation at break and the percent water sorption were increased with the increasing of the percentage of dibutyl phthalate as plasticizer whereas the tensile strength were decreased. The release rate of the drug from capsules containing granules decreased when increasing the amount of aqueous polymeric dispersion. The type of aqueous polymeric dispersion affected the release of drug. Capsules containing larger granules exhibited slower release of drug than smaller granules. For similar of drug release, the amount of Eudragit[®] NE 30D used in coating was higher than that of the Surelease[®]. The drug release model from capsules containing granules fitted mostly to the first-order model. Some formulations of coated granules exhibited drug release conformed to the USP.

ภาควิชา เภสัชศาสตร์
สาขาวิชา เภสัชกรรมคลินิก
ปีการศึกษา ๒๕๓๙

ลายมือชื่อนิสิต ลักษณ์ ตั้งกานดา
ลายมือชื่ออาจารย์ที่ปรึกษา Darnpimol C. Ritthidej
ลายมือชื่ออาจารย์ที่ปรึกษาร่วม Wichein Thanindratarn

Acknowledgments

I would like to express my sincere gratitude to my thesis advisor and coadvisor, Associate Professor Gampimol C. Rithidej, Ph.D. and Assistant Professor Wichein Thanindrataarn, M.Sc. in Pharm, for the valuable advice, guidance and encouragement throughout this study. The patience, kindness and understanding are also deeply appreciated.

A special appreciation is also given to the Graduate School, Chulalongkorn University for granting partial financial support to fulfill this investigation.

To the other members of thesis committee, I wish to appreciate for their valuable suggestions.

The special acknowledgement is given to all staffs in the Department of Manufacturing Pharmacy for their assistance.

Finally, I would like to express my infinite thanks and deepest gratitude to my parents for their endless love, continuous support, care and understanding.

CONTENTS

	Page
Thai Abstract.....	iv
English Abstract.....	v
Acknowledgement.....	vi
List of Figures.....	viii
List of Tables.....	xvi
 Chapter	
I Introduction.....	1
II Experiment.....	37
III Results.....	54
IV Discussion and Conclusions.....	143
References.....	159
Appendix.....	164
Vitae.....	232

สถาบันวิทยบริการ
จุฬาลงกรณ์มหาวิทยาลัย

LIST OF FIGURES

Figure		Page
1	Chemical structure of theophylline.....	10
2	Structural formula of Surelease® (ethylcellulose aqueous dispersion).....	11
3	Structural formula of methacrylate ester copolymers.....	16
4	Top spray coater.....	22
5	Wurster bottom spray coater.....	24
6	Rotor tangential spary coater.....	26
7	The photograph of fluidized bed coater.....	42
8	The potograph of standard template of aqueous polymeric films.....	43
9	The photomicrographs of uncoated theophylline granules of 16/18 mesh size.....	56
10	The photomicrographs of uncoated theophylline granules of 18/20 mesh size.....	57
11	The photomicrographs of uncoated theophylline granules of 20/25 mesh size.....	58
12	The photomicrographs of uncoated theophylline granules of 25/30 mesh size.....	59
13	The photomicrographs of 1.92% Surelease® coated granules of 16/18 mesh size.....	61
14	The photomicrographs of 3.05% Surelease® coated granules of 16/18 mesh size.....	62
15	The photomicrographs of 3.87% Surelease® coated granules of 16/18 mesh size.....	63

Figure	Page
16 The photomicrographs of 6.29% Surelease® coated granules of 16/18 mesh size.....	64
17 The photomicrographs of 2.04% Surelease® coated granules of 18/20 mesh size.....	65
18 The photomicrographs of 3.29% Surelease® coated granules of 18/20 mesh size.....	66
19 The photomicrographs of 7.61% Surelease® coated granules of 18/20 mesh size.....	67
20 The photomicrographs of 9.00% Surelease® coated granules of 18/20 mesh size.....	68
21 The photomicrographs of 2.48% Surelease® coated granules of 20/25 mesh size.....	69
22 The photomicrographs of 3.86% Surelease® coated granules of 20/25 mesh size.....	70
23 The photomicrographs of 8.39% Surelease® coated granules of 20/25 mesh size.....	71
24 The photomicrographs of 12.12% Surelease® coated granules of 20/25 mesh size.....	72
25 The photomicrographs of 7.84% Surelease® coated granules of 25/30 mesh size.....	73
26 The photomicrographs of 12.18% Surelease® coated granules of 25/30 mesh size.....	74
27 The photomicrographs of 13.40% Surelease® coated granules of 25/30 mesh size.....	75
28 The photomicrographs of 17.05% Surelease® coated granules of 25/30 mesh size.....	76

Figure		Page
29	The photomicrographs of 2.51% Eudragit® NE 30D coated granules of 16/18 mesh size.....	77
30	The photomicrographs of 5.96% Eudragit® NE 30D coated granules of 16/18 mesh size.....	78
31	The photomicrographs of 8.53% Eudragit® NE 30D coated granules of 16/18 mesh size.....	79
32	The photomicrographs of 14.76% Eudragit® NE 30D coated granules of 16/18 mesh size.....	80
33	The photomicrographs of 4.04% Eudragit® NE 30D coated granules of 18/20 mesh size.....	81
34	The photomicrographs of 6.75% Eudragit® NE 30D coated granules of 18/20 mesh size.....	82
35	The photomicrographs of 10.93% Eudragit® NE 30D coated granules of 18/20 mesh size.....	83
36	The photomicrographs of 19.12% Eudragit® NE 30D coated granules of 18/20 mesh size.....	84
37	The photomicrographs of 5.56% Eudragit® NE 30D coated granules of 20/25 mesh size.....	85
38	The photomicrographs of 11.46% Eudragit® NE 30D coated granules of 20/25 mesh size.....	86
39	The photomicrographs of 15.00% Eudragit® NE 30D coated granules of 20/25 mesh size.....	87
40	The photomicrographs of 20.76% Eudragit® NE 30D coated granules of 20/25 mesh size.....	88
41	The photomicrographs of 8.83% Eudragit® NE 30D coated granules of 25/30 mesh size.....	89

Figure		Page
42	The photomicrographs of 14.57% Eudragit® NE 30D coated granules of 25/30 mesh size.....	90
43	The photomicrographs of 17.68% Eudragit® NE 30D coated granules of 25/30 mesh size.....	91
44	The photomicrographs of 22.37% Eudragit® NE 30D coated granules of 25/30 mesh size.....	92
45	Influence of dibutyl phthalate concentration on percent elongation at break of Surelease® films.....	105
46	Influence of dibutyl phthalate concentration on tensile strength of Surelease® films.....	106
47	Influence of dibutyl phthalate concentration on percentage of water sorption of Surelease® films.....	108
48	A. The release profiles of uncoated theophylline granules of various sizes.....	108
	B. The release rate profiles of uncoated theophylline granules of various sizes.....	110
49	A. The release profiles of Surelease® coated granules of 16/18 mesh sizes.....	110
	B. The release rate profiles of Surelease® coated granules of 16/18 mesh sizes.....	110
50	A. The release profiles of Surelease® coated granules of 18/20 mesh sizes.....	111
	B. The release rate profiles of Surelease® coated granules of 18/20 mesh sizes.....	111

Figure		Page
51	A. The release profiles of Surelease® coated granules of 20/25 mesh sizes.....	112
	B. The release rate profiles of Surelease® coated granules of 20/25 mesh sizes.....	112
52	A. The release profiles of Surelease® coated granules of 25/30 mesh sizes.....	113
	B. The release rate profiles of Surelease® coated granules of 25/30 mesh sizes.....	113
53	A. The release profiles of Eudragit® NE 30D coated granules of 16/18 mesh sizes.....	115
	B. The release rate profiles of Eudragit® NE 30D coated granules of 16/18 mesh sizes.....	115
54	A. The release profiles of Eudragit® NE 30D coated granules of 18/20 mesh sizes.....	116
	B. The release rate profiles of Eudragit® NE 30D coated granules of 18/20 mesh sizes.....	116
55	A. The release profiles of Eudragit® NE 30D coated granules of 20/25 mesh sizes.....	117
	B. The release rate profiles of Eudragit® NE 30D coated granules of 20/25 mesh sizes.....	117
56	A. The release profiles of Eudragit® NE 30D coated granules of 25/30 mesh sizes.....	118
	B. The release rate profiles of Eudragit® NE 30D coated granules of 25/30 mesh sizes.....	118
57	A. The release profiles of Theo-24®.....	119
	B. The release rate profiles of Theo-24®.....	119

Figure	Page
58 Influence of various sizes (18/20 and 20/25 mesh) in about 2% coating level of Surelease® on the theophylline release profiles from coated granules.....	120
59 Influence of various sizes (16/18, 18/20 and 20/25 mesh) in about 3% coating level of Surelease® on the theophylline release profiles from coated granules.....	120
60 Influence of various sizes (18/20 and 25/30 mesh) in about 7% coating level of Surelease® on the theophylline release profiles from coated granules.....	121
61 Influence of various sizes (20/25 and 25/30 mesh) in about 12% coating level of Surelease® on the theophylline release profiles from coated granules.....	121
62 Influence of various sizes (16/18 and 20/25 mesh) in about 5% coating level of Eudragit® NE 30D on the theophylline release profiles from coated granules.....	123
63 Influence of various sizes (16/18 and 25/30 mesh) in about 8% coating level of Eudragit® NE 30D on the theophylline release profiles from coated granules.....	123
64 Influence of various sizes (16/18 and 25/30 mesh) in about 14% coating level of Eudragit® NE 30D on the theophylline release profiles from coated granules.....	124
65 The release profiles between 6.29% Surelease® coated granules and 14.76% Eudragit® NE 30D coated granules for granules of 16/18 mesh.....	126
66 The release profiles between 7.61% Surelease® coated granules and 19.12% Eudragit® NE 30D coated granules for granules of 18/20 mesh.....	126

Figure	Page
67 The release profiles between 3.86% Surelease® coated granules and 5.56% Eudragit® NE 30D coated granules for granules of 20/25 mesh.....	127
68 The release profiles between 7.84% Surelease® coated granules and 14.57% Eudragit® NE 30D coated granules for granules of 25/30 mesh.....	127
69 The release profiles from all selected formulations of Surelease® coated granules compared with Theo-24®	129
70 The release profiles from 6.29% Surelease® coated granules of 16/18 mesh compared with Theo-24®	130
71 The release profiles from 7.61% Surelease® coated granules of 18/20 mesh compared with Theo-24®	130
72 The release profiles from 8.39% Surelease® coated granules of 20/25 mesh compared with Theo-24®	131
73 The release profiles from 12.18% and 13.40% Surelease® coated granules of 25/30 mesh compared with Theo-24®	131
74 The release profiles from all selected formulations of Eudragit® NE 30D coated granules compared with Theo-24®	132
75 The release profiles from 14.76% Eudragit® NE 30D coated granules of 16/18 mesh compared with Theo-24®	133
76 The release profiles from 19.12% Eudragit® NE 30D coated granules of 18/20 mesh compared with Theo-24®	133
77 The release profiles from 20.76% Eudragit® NE 30D coated granules of 20/25 mesh compared with Theo-24®	134
78 The release profiles from 22.37% Eudragit® NE 30D coated granules of 25/30 mesh compared with Theo-24®	134

Figure		Page
79	Diffusion control of drug release through a aqueous polymeric film.....	151
80	Structural formula of Surelease®.....	153
81	Structural formula of Eudragit® NE 30D.....	154
82	Calibration curve of theophylline in phosphate buffer pH 6.6	166
83	Calibration curve of HPLC method.....	166
84	The representative chromatogram of the construction of calibration curve.....	167
	A. Conc. of theophylline granules = 20.0 $\mu\text{g/ml}$	
	B. Conc. of theophylline granules = 140.0 $\mu\text{g/ml}$	



**สถาบันวิทยบริการ
จุฬาลงกรณ์มหาวิทยาลัย**

LIST OF TABLES

Table		Page
1	Characteristics of three fluidized bed coating process.....	27
2	Composition of theophylline granules for coating.....	40
3	Composition of ethylcellulose aqueous polymeric coating.....	42
4	Composition of poly(ethylacrylate methylmethacrylate) aqueous polymeric coating.....	44
5	The coating conditons by top spray method.....	44
6	The amount of Surelease® for coating a batch of theophylline granules (350 g).....	44
7	The amount of Eudragit® NE 30D for coating a batch of theophylline granules (350 g).....	45
8	The percent coating levels of Surelease® on granules of different sizes.....	54
9	The percent coating levels of Eudragit® NE 30D on granules of different sizes.....	55
10	Physical properties of uncoated theophylline granules.....	98
11	Mean sizes of Surelease® coated granules.....	98
12	Mean sizes of Eudragit® NE 30D coated granules.....	99
13	Bulk densities of Surelease® coated granules.....	99
14	Bulk densities of Eudragit® NE 30D coated granules.....	99
15	Tapped densities of Surelease® coated granules.....	100
16	Tapped densities of Eudragit® NE 30D coated granules.....	100
17	Carr's compressibilities of Surelease® coated granules.....	100
18	Carr's compressibilities of Eudragit® NE 30D coated granules...	101
19	Flow rates of Surelease® coated granules.....	101

Table		Page
20	Flow rates of Eudragit® NE 30D coated granules.....	101
21	Angles of repose of Surelease® coated granules.....	102
22	Angles of repose of Eudragit® NE 30D coated granules.....	102
23	Moisture contents of Surelease® coated granules.....	102
24	Moisture contents of Eudragit® NE 30D coated granules.....	103
25	The angle of repose of the granules.....	103
26	Specific surface area of uncoated granules.....	103
27	Specific surface area of Surelease® coated granules.....	104
28	Specific surface area of Eudragit® NE 30D coated granules.....	104
29	Percent elongation at break and tensile strength of aqueous polymeric films.....	104
30	Percent water sorption of aqueous polymeric films.....	107
31	USP XXIII requirement and cumulative percent release of drug from selected formulations of Surelease® coated granules.....	137
32	USP XXIII requirement and cumulative percent release of drug from selected formulations of Eudragit® NE 30D coated granules.....	138
33	Correlation coefficient of the relationships between percent drug released versus time (A), percent drug released versus square root time (B), and log percent drug remained versus time (C) from uncoated granules.....	138
34	Correlation coefficient of the relationships between percent drug released versus time (A), percent drug released versus square root time (B), and log percent drug remained versus time (C) from theophylline granules coating with Surelease®.....	139

Table		Page
35	Correlation coefficient of the relationships between percent drug released versus time (A), percent drug released versus square root time (B), and log percent drug remained versus time (C) from theophylline granules coating with Eudragit® NE 30D	140
36	Correlation coefficient of the relationships between percent drug released versus time (A), percent drug released versus square root time (B), and log percent drug remained versus time (C) from theophylline granules coating with Theo-24® (Commercial Product).....	141
37	Comparison of linearity between plots of rate of release against reciprocal amount (1/Q and amount (Q) of theophylline release from theophylline granules coating with Surelease®	141
38	Comparison of linearity between plots of rate of release against reciprocal amount (1/Q and amount (Q) of theophylline release from theophylline granules coating with Eudragit® NE 30D.....	142
39	Comparison of linearity between plots of rate of release against reciprocal amount (1/Q and amount (Q) of theophylline release from Theo-24®, 200 mg.....	142
40	Absorbance of theophylline in phosphate buffer pH 6.6 determine at 272 nm.....	164
41	Ratio between peak area of internal standard and theophylline anhydrous with HPLC method at 272 nm.....	165
42	Amount percent of theophylline release from uncoated granules in phosphate buffer pH 6.6.....	168

Table		Page
43	Amount percent of theophylline release from 16/18 mesh size of theophylline granules coating with Surelease® in phosphate buffer pH 6.6.....	169
44	Amount percent of theophylline release from 18/20 mesh size of theophylline granules coating with Surelease® in phosphate buffer pH 6.6.....	170
45	Amount percent of theophylline release from 20/25 mesh size of theophylline granules coating with Surelease® in phosphate buffer pH 6.6.....	171
46	Amount percent of theophylline release from 25/30 mesh size of theophylline granules coating with Surelease® in phosphate buffer pH 6.6.....	172
47	Amount percent of theophylline release from 16/18 mesh size of theophylline granules coating with Eudragit® NE 30D in phosphate buffer pH 6.6.....	173
48	Amount percent of theophylline release from 18/20 mesh size of theophylline granules coating with Eudragit® NE 30D in phosphate buffer pH 6.6.....	174
49	Amount percent of theophylline release from 20/25 mesh size of theophylline granules coating with Eudragit® NE 30D in phosphate buffer pH 6.6.....	175
50	Amount percent of theophylline release from 25/30 mesh size of theophylline granules coating with Eudragit® NE 30D in phosphate buffer pH 6.6.....	176
51	Amount percent of theophylline release from Theo-24®, 200 mg, in phosphate buffer pH 6.6.....	177

Table		Page
52	The release rate of uncoated granules in phosphate buffer pH 6.6.....	178
53	The release rate of 16/18 mesh size of theophylline granules coating with Surelease® in phosphate buffer pH 6.6.....	179
54	The release rate of 18/20 mesh size of theophylline granules coating with Surelease® in phosphate buffer pH 6.6.....	180
55	The release rate of 20/25 mesh size of theophylline granules coating with Surelease® in phosphate buffer pH 6.6.....	181
56	The release rate of 25/30 mesh size of theophylline granules coating with Surelease® in phosphate buffer pH 6.6.....	182
57	The release rate of 16/18 mesh size of theophylline granules coating with Eudragit® NE 30D in phosphate buffer pH 6.6.....	183
58	The release rate of 18/20 mesh size of theophylline granules coating with Eudragit® NE 30D in phosphate buffer pH 6.6.....	184
59	The release rate of 20/25 mesh size of theophylline granules coating with Eudragit® NE 30D in phosphate buffer pH 6.6.....	185
60	The release rate of 25/30 mesh size of theophylline granules coating with Eudragit® NE 30D in phosphate buffer pH 6.6.....	186
61	The release rate of commercial product, Theo-24® 200 mg, in phosphate buffer pH 6.6.....	187
62	Values of rate (dQ/dt), amount released (Q), and the corresponding reciprocal (1/Q) for the release of 16/18 mesh size of theophylline granules coating with Surelease®	188
63	Values of rate (dQ/dt), amount released (Q), and the corresponding reciprocal (1/Q) for the release of 18/20 mesh size of theophylline granules coating with Surelease®	189

Table	Page
64 Values of rate (dQ/dt), amount released (Q), and the corresponding reciprocal ($1/Q$) for the release of 20/25 mesh size of theophylline granules coating with Surelease®.....	190
65 Values of rate (dQ/dt), amount released (Q), and the corresponding reciprocal ($1/Q$) for the release of 25/30 mesh size of theophylline granules coating with Surelease®.....	191
66 Values of rate (dQ/dt), amount released (Q), and the corresponding reciprocal ($1/Q$) for the release of 16/18 mesh size of theophylline granules coating with Eudragit® NE 30D....	192
67 Values of rate (dQ/dt), amount released (Q), and the corresponding reciprocal ($1/Q$) for the release of 18/20 mesh size of theophylline granules coating with Eudragit® NE 30D....	193
68 Values of rate (dQ/dt), amount released (Q), and the corresponding reciprocal ($1/Q$) for the release of 20/25 mesh size of theophylline granules coating with Eudragit® NE 30D....	194
69 Values of rate (dQ/dt), amount released (Q), and the corresponding reciprocal ($1/Q$) for the release of 25/30 mesh size of theophylline granules coating with Eudragit® NE 30D....	195
70 Values of rate (dQ/dt), amount released (Q), and the corresponding reciprocal ($1/Q$) for the release of commercial product, Theo-24® 200 mg.....	196
71 comparison of percent water sorption between Surelease® Films and Eudragit® NE 30D film.....	196
72 The t-values of percent water sorption.....	197
73 The t-values of percentage drug release for various levels of Surelease®	197

Table		Page
74	The t-values of percentage drug release between Surelease® and Eudragit® NE 30D coated granules	198
75	The t-values of percentage drug release between Surelease® coated granules and commercial product.....	199
76	Comparison of percent compressibility among various sizes of uncoated theophylline granules.....	200
77	The F-ratio of percent compressibility for uncoated theophylline granules of various sizes (data from Table 76).....	200
78	Comparison of flow rates among various sizes of uncoated theophylline granules.....	201
79	The F-ratio of flow rates for uncoated theophylline granules of various sizes (data from Table 78).....	201
80	Comparison of percent compressibility among various percent coating levels of Surelease® for granules of 16/18 mesh size.....	202
81	Comparison of percent compressibility among various percent coating levels of Surelease® for granules of 18/20 mesh size.....	202
82	Comparison of percent compressibility among various percent coating levels of Surelease® for granules of 20/25 mesh size.....	203
83	Comparison of percent compressibility among various percent coating levels of Surelease® for granules of 25/30 mesh size.....	203
84	Comparison of percent compressibility among various percent coating levels of Eudragit® NE 30D for granules of 16/18 mesh size.....	204
85	Comparison of percent compressibility among various percent coating levels of Eudragit® NE 30D for granules of 18/20 mesh size.....	204

Table		Page
86	Comparison of percent compressibility among various percent coating levels of Eudragit® NE 30D for granules of 20/25 mesh size.....	205
87	Comparison of percent compressibility among various percent coating levels of Eudragit® NE 30D for granules of 25/30 mesh size.....	205
88	Comparison of flow rates among various percent coating levels of Surelease® for granules of 16/18 mesh size.....	206
89	Comparison of flow rates among various percent coating levels of Surelease® for granules of 18/20 mesh size.....	206
90	Comparison of flow rates among various percent coating levels of Surelease® for granules of 20/25 mesh size.....	207
91	Comparison of flow rates among various percent coating levels of Surelease® for granules of 25/30 mesh size.....	207
92	Comparison of flow rates among various percent coating levels of Eudragit® NE 30D for granules of 16/18 mesh size.....	208
93	Comparison of flow rates among various percent coating levels of Eudragit® NE 30D for granules of 18/20 mesh size.....	208
94	Comparison of flow rates among various percent coating levels of Eudragit® NE 30D for granules of 20/25 mesh size.....	209
95	Comparison of flow rates among various percent coating levels of Eudragit® NE 30D for granules of 25/30 mesh size.....	209
96	The F-ratio of percent compressibility for Surelease® coated granules of various sizes (data from Tables 80-83).....	210
97	The F-ratio of percent compressibility for Eudragit® NE 30D coated granules of various sizes (data from Tables 84-87).....	211

Table		Page
98	The F-ratio of flow rates for Surelease® coated granules of various sizes (data from Tables 88-91).....	212
99	The F-ratio of flow rates for Eudragit® NE 30D coated granules of various sizes (data from Tables 92-95).....	213
100	The t-values of percent compressibility for uncoated theophylline granules of various sizes (data from Tables 76).....	214
101	The t-values of flow rates for uncoated theophylline granules of various sizes (data from Table 78).....	214
102	The t-values of percent compressibility for surelease® coated granules of 16/18 mesh size (data from Table 80).....	215
103	The t-values of percent compressibility for Surelease® coated granules of 18/20 mesh size (data from Table 81).....	216
104	The t-values of percent compressibility for Eudragit® NE 30D coated granules of 16/18 mesh size (data from Table 84).....	217
105	The F-ratio of percent compressibility between uncoated and Surelease® coated granules of various sizes (data from Tables 76 and 80-83).....	218
106	The F-ratio of percent compressibility between uncoated and Eudragit® NE 30D coated granules of various sizes (data from Tables 76 and 84-87).....	219
107	The F-ratio of flow rates between uncoated and Surelease® coated granules of various sizes (data from Tables 78 and 88-91).....	220
108	The F-ratio of flow rates between uncoated and Eudragit® NE 30D coated granules of various sizes (data from Tables 78 and 92-95).....	221

Table		Page
109	The t-values of percent compressibility between uncoated and Surelease® coated granules of 16/18 mesh size (data from Tables 76 and 80).....	222
110	The t-values of percent compressibility between uncoated and Surelease® coated granules of 18/20 mesh size (data from Tables 76 and 81).....	223
111	The t-values of percent compressibility between uncoated and Eudragit® NE 30D coated granules of 16/18 mesh size (data from Tables 76 and 84).....	224
112	The t-values of percent compressibility between uncoated and Eudragit® NE 30D coated granules of 20/25 mesh size (data from Tables 76 and 86).....	225
113	The t-values of percent compressibility between uncoated and Eudragit® NE 30D coated granules of 25/30 mesh size (data from Tables 76 and 87).....	226
114	The t-values of flow rates between uncoated and Surelease® coated granules of 18/20 mesh size (data from Tables 78 and 89).....	227
115	The t-values of flow rates between uncoated and Surelease® coated granules of 20/25 mesh size (data from Tables 78 and 89).....	228
116	The t-values of flow rates between uncoated and Eudragit® NE 30D coated granules of 18/20 mesh size (data from Tables 78 and 93).....	229
117	The t-values of flow rates between uncoated and Eudragit® NE 30D coated granules of 20/25 mesh size (data from Tables 78 and 94).....	230

Table		Page
118	The t-values of flow rates between uncoated and Eudragit® NE 30D coated granules of 25/30 mesh size (data from Tables 78 and 95).....	231



ABBREVIATIONS

bar	N/cm ²
°C	degree celcius
cm	centimeter
Conc.	Concentration
g	gram
HCl	Hydrochloric acid
hr.	hour
mg	milligram
min.	minute
ml	milliliter
ml/min	milliliter per minute
N	Normality
NaOH	Sodium hydroxide
nm	nanometer
SD	Standard deviation
sec.	second
SEM	Scanning Electron Microscope
UV	ultraviolet
wt.	weight
µg/ml	microgram per milliliter
µm (µ)	micrometer